

AccuVert™ HCV Seroconversion Panel

0610-0223 / Batch #10683758


OVERVIEW

AccuVert™ HCV Seroconversion Panel 0610-0223 / Batch #10683758 is a 6-member validation panel of undiluted, naturally occurring plasma samples (1 vial per member, 1.2 mL per vial). Panel members represent serial bleeds collected over the course of 88 days in 2008 from a 46-year old male with HCV genotype 2b. Each sample represents a single collection event. No preservatives were added.

Test results from commercially available HCV assays are included for characterization of the panel members. This panel of human plasma samples converts from negative to positive for HCV Antibody.

For Research Use Only. Not for use in diagnostic procedures. Data are offered for informational purposes. LGC SeraCare Life Sciences does not claim that others can duplicate test results exactly.

CAUTION: Potentially infectious materials. Follow Universal Precautions. Some panel members were found positive for anti-HCV; all were found negative for anti-HIV-1/2 and HBsAg. This does not ensure the absence of these or other human pathogens.



This graph demonstrates the evolution of early HCV infection amongst panel members from the Abbott ARCHITECT (Anti-HCV) and the Roche COBAS (HCV RNA) test methods.

AccuVert™ HCV Seroconversion Panel

0610-0223 / Batch #10683758

Panel Member Information

Panel Member	SeraCare Batch #	SeraCare Donor ID	Bleed Date	Days Since 1 st Bleed
01	9235383	BD107287	22-Sep-2008	0
02	9235384	BD107287	25-Sep-2008	3
03	9235385	BD107287	29-Sep-2008	7
04	9235300	BD107287	20-Nov-2008	59
05	9235519	BD107287	25-Nov-2008	64
06	9237286	BD107287	19-Dec-2008	88

HCV RNA¹

Panel Member	Roche cobas® 6800/8800 HCV RNA (IU/mL)	Abbott Alinity m HCV RNA (IU/mL)
01	1.21 x 10⁶	9.80 x 10⁵
02	4.74 x 10⁵	3.79 x 10⁵
03	4.10 x 10⁵	2.98 x 10⁵
04	3.23 x 10³	2.84 x 10³
05	1.26 x 10³	1.20 x 10³
06	6.30 x 10²	4.16 x 10²
Test Date	12-Oct-2023	12-Oct-2023
Test Site	RL	RL
Kit Part Code	NA	NA
Kit Lot No.	NA	385356
Kit Exp. Date	NA	10-Feb-2025

¹Results are reported as the mean result of duplicate testing; positive/reactive results are noted in bold red.
RL = Reference Lab; NA = Not Available

AccuVert™ HCV Seroconversion Panel

0610-0223 / Batch #10683758

Anti-HCV¹

Panel Member	Abbott ARCHITECT Anti-HCV (S/CO)	Roche cobas® Elecsys Anti-HCV II (COI)	Abbott Alinity s Anti-HCV (S/CO)	ORTHO HCV Version 3.0 ELISA (S/CO)
01	0.40	0.08	0.20	0.05
02	0.35	0.08	0.18	0.04
03	0.27	0.09	0.14	0.03
04	9.24	176.50	5.31	4.71
05	10.89	182.35	5.92	>6.56
06	13.37	87.54	10.18	>6.56
Test Date	10-Oct-2023	08-Nov-2023	10-Oct-2023	10-Oct-2023
Test Site	LGC	RL	RL	LGC
Kit Part Code	1L79	NA	NA	6195566
Kit Lot No.	52011BE00	NA	51608BE00	SXE037
Kit Exp. Date	05-Mar-2024	NA	04-Apr-2024	20-Jan-2024

¹Results are reported as the mean result of duplicate testing; positive/reactive results are noted in bold red.
RL = Reference Lab; LGC = LGC SeraCare; NA = Not Available

HCV Confirmatory Assay¹

INNOGENETICS® INNO-LIA™ HCV Score

Band Results

Panel Member	HCV Antibody Interpretation	C1	C2	E2	NS3	NS4	NS5
01	NEG	-	-	-	-	-	-
02	NEG	-	-	-	-	-	-
03	NEG	-	-	-	-	-	-
04	POS	1+	-	-	2+	-	+/-
05	POS	2+	-	-	2+	-	+/-
06	POS	2+	-	-	3+	-	-
Test Date	20-Oct-2023						
Test Site	RL						
Kit Part Code	NA						
Kit Lot No.	411025						
Kit Exp. Date	30-Jun-2024						

¹Results are reported as singlicate testing; positive/reactive results are noted in bold red.
RL = Reference Lab; NA = Not Available
NEG = Negative; POS = Positive

The package insert for this panel can be found at
www.seracare.com.

A printed copy of the package insert or data sheet may be requested by email at CDx-Info@LGCGroup.com or by phone at 508.244.6400.