

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified PRB109(M) (0800-0301) / Batch #10311491

OVERVIEW


AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified PRB109(M) (0800-0301) / Batch #10311491 is a 15-member validation panel of undiluted, naturally occurring plasma samples (1 vial per member, 1.0 mL per vial). Panel members represent bleeds from multiple individuals positive for antigen/antibodies to HIV. Each sample represents a single collection event. No preservatives were added.

Test results from commercially-available HIV assays are included for characterization of the panel members. This panel of human plasma samples demonstrates a range of antibody reactivity which approaches the sensitivity limits for several HIV test methods. One sample is included as a non-reactive sample and is negative for all HIV test methods performed.

For Research Use Only. Not for use in diagnostic procedures. Data are provided for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.

CAUTION: Potentially infectious materials. Follow Universal Precautions. The units that make up this panel were tested and found negative for anti-HCV and HBsAg. This does not ensure the absence of these or other human pathogens.

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified



This graph demonstrates HIV-1 Ag/Ab reactivity amongst panel members utilizing test results from the Abbott Prism HIV Ag/Ab Combo method. Panel members 5, 6, 9, 18 and 20 are no longer available. Data are provided for informational purposes only.

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified PRB109(M) (0800-0301) / Batch #10311491

Panel Member Information

Panel Member	SeraCare Batch #	SeraCare Donor ID #	Bleed Date	Fiebig Stage ¹
01	BM202249	NA	06-Nov-01	IV
02	BM202261	NA	18-Jan-02	II / III
03	BM202265	NA	18-Feb-02	IV
04	BM202266	NA	20-Feb-02	IV
05	BM202267	NA	02-May-02	IV
06	BM202271	NA	06-Jun-02	IV
07	BM202272	NA	31-May-02	III
08	9178194	NA	31-Jul-03	NA
09	BM202291	NA	22-Oct-02	IV
10	BM206239	BD102670	22-May-06	IV
11	BM206240	BD102670	24-May-06	IV
12	BM202161	BD100797	21-Oct-04	V
13	BM217811	BD105709	14-Sep-07	IV
14	BM217812	BD105709	19-Sep-07	IV
15	BM217813	BD105709	21-Sep-07	V
16	9222531	BD105878	20-Dec-07	IV
17	BM202297	NA	12-Oct-02	IV
18	BM202296	NA	14-Oct-02	V
19	9230846	BD106689	29-Apr-08	V
20	BM202278	NA	13-Jun-02	V

¹Fiebig stage is assigned according to EW Fiebig et al, Dynamics of HIV viremia and antibody seroconversion in plasma donors: implications for diagnosis and staging of primary HIV infection.

Panel members 5, 6, 9, 18 and 20 are no longer available. Data are provided for informational purposes only.

NA = Not Available

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified PRB109(M) (0800-0301) / Batch #10311491

HIV RNA and Antigen

Panel Member	Roche Amplicor Monitor HIV-1 RNA (c/mL)	Perkin Elmer HIV Antigen (s/co) ¹	Zeptometrix HIV Antigen (s/co) ¹
01	5.5 x 10⁵	33.2	6.2
02	4.8 x 10⁶	>46.0	24.2
03	2.3 x 10⁵	5.3	1.5
04	2.7 x 10⁵	10.5	2.5
05	3.2 x 10⁴	1.4	1.0
06	2.3 x 10⁵	12.6	3.6
07	6.2 x 10⁴	3.6	1.3
08	BLD	0.3	0.7
09	2.1 x 10⁵	5.3	2.4
10	1.2 x 10⁵	3.5	1.3
11	5.7 x 10⁵	10.6	2.4
12	4.1 x 10⁵	8.7	2.4
13	5.0 x 10⁵	19.7	4.2
14	2.4 x 10⁶	>46.0	11.8
15	7.0 x 10⁵	>46.0	8.4
16	5.7 x 10⁶	>46.0	16.5
17	1.7 x 10⁵	28.0	6.8
18	7.8 x 10⁴	3.7	1.3
19	2.1 x 10⁴	0.4	0.7
20	6.9 x 10³	1.1	0.9
Test Date	18-Jun-09, 29-Jun-09	18-Jun-09	16-Jun-09
Test Site	SC	SC	SC
Kit Part Code	83377	NEK 050B	0801111
Kit Lot No.	K10107	990-090302	304824
Kit Exp. Date	31-Oct-09	01-Jan-10	15-Dec-09
Kit Regulatory Status	IVD/CE	RUO	RUO

¹Ratios ≥1.0 are considered positive/reactive and noted in bold red.

Panel members 5, 6, 9, 18 and 20 are no longer available. Data are provided for informational purposes only.

SC = SeraCare

IVD = In Vitro Diagnostic; CE = Conformité Européenne or CE Marking; RUO = Research Use Only

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified

PRB109(M) (0800-0301) / Batch #10311491

HIV Antibody (s/co)^{1,2}

Panel Member	Abbott HIVAB HIV-1/HIV-2 (rDNA)	Bio-Rad GS HIV-1/HIV-2 Plus O	Genetic Systems HIV-1 rLAV
	EIA	EIA	EIA
01	4.1	>10.9	0.2
02	0.7	>10.9	0.1
03	0.6	8.5	0.1
04	2.1	>10.9	0.2
05	18.4	>10.9	1.7
06	3.6	4.1	1.2
07	0.7	1.3	0.2
08	0.1	0.1	0.1
09	2.5	>10.9	0.6
10	10.3	>10.9	1.0
11	7.4	>10.9	2.1
12	11.8	>10.9	2.1
13	0.5	0.3	0.2
14	10.1	>10.9	3.1
15	7.2	>10.9	3.9
16	6.1	>10.9	0.3
17	2.2	>10.9	0.4
18	7.7	>10.9	1.5
19	16.3	>10.9	3.3
20	1.5	>10.9	2.2
Test Date	15 Jun 09	15 Jun 09	15 Jun 09
Test Site	SC	SC	SC
Kit Part Code	3A77	32588	32511
Kit Lot No.	76625M100	324CBB	221CAA
Kit Exp. Date	02 Aug 09	13 Aug 09	01 Jul 09
Kit Regulatory Status	IVD	IVD	RUO

¹EIA results are reported as a signal to cutoff ratio (s/co); Ratios ≥ 1.0 are considered positive/reactive and noted in bold red.

²Results are reported as the mean result of duplicate testing.

Panel members 5, 6, 9, 18 and 20 are no longer available. Data are provided for informational purposes only.

SC = SeraCare

IVD = In Vitro Diagnostic; RUO = Research Use Only

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified PRB109(M) (0800-0301) / Batch #10311491

HIV Antigen/Antibody (s/co)¹

Panel Member	Abbott ARCHITECT Combo HIV Ag/Ab	Abbott PRISM Combo HIV Ag/Ab	Murex Combination EIA HIV Ag/Ab	Roche Elecsys Combi HIV Ag/Ab
01	40.9	16.6	17.8	8.2
02	281.3	121.7	18.9	36.7
03	3.5	2.3	2.8	5.7
04	10.7	4.8	6.4	2.0
05	10.2	4.9	17.5	37.5
06	4.1	4.1	3.8	0.9
07	3.9	2.7	3.4	1.1
08	0.3	0.2	0.4	0.3
09	10.6	4.6	8.0	4.1
10	20.5	8.4	12.3	9.5
11	36.5	15.5	14.3	50.1
12	34.1	7.6	17.5	193.0
13	23.6	12.6	9.6	1.8
14	127.5	50.3	18.8	41.9
15	127.9	41.9	17.8	108.0
16	116.7	55.1	17.9	24.0
17	41.6	16.9	15.6	5.6
18	58.6	6.7	14.3	13.2
19	23.5	6.0	17.9	291.3
20	13.7	3.4	7.9	39.9
Test Date	29-Jul-09	05-Aug-09	03-Aug-09, 07-Aug-09	13-Aug-09
Test Site	RL	RL	RL	RL
Kit Part Code	4J27-20	7G46-48	GE42 (7G79-02)	4860446
Kit Lot No.	72331HN00	79496HN00	L283810	153019-01
Kit Exp. Date	14-Sep-09	28-Sep-09	31-Mar-10	31-Oct-09
Kit Regulatory Status	IVD/CE	IVD/CE	IVD/CE	IVD/CE

¹Results are reported as a signal to cutoff ratio (s/co); Ratios ≥ 1.0 are considered positive/reactive and noted in bold red.

Panel members 5, 6, 9, 18 and 20 are no longer available. Data are provided for informational purposes only.

RL = Reference Lab

IVD = In Vitro Diagnostic; CE = Conformité Européenne or CE Marking

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified

PRB109(M) (0800-0301) / Batch #10311491

HIV Rapid Tests (s/co)¹

Panel Member	Bio-Rad Multispot HIV-1/HIV-2		OraSure OraQuick ADVANCE® HIV-1/2 Result	Inverness Determine® HIV-1/2 Ag/Ab Combo	
	HIV-1 Result	HIV-2 Result		HIV Ag Result	HIV-1/2 Ab Result
01	NEG	NEG	NEG	POS	NEG
02	NEG	NEG	NEG	POS	POS
03	NEG	NEG	NEG	NEG	NEG
04	NEG	NEG	NEG	NEG	POS
05	NEG	NEG	NEG	NEG	POS
06	NEG	NEG	NEG	POS	NEG
07	NEG	NEG	NEG	NEG	NEG
08	NEG	NEG	NEG	NEG	NEG
09	NEG	NEG	NEG	NEG	POS
10	NEG	NEG	NEG	NEG	POS
11	NEG	NEG	POS	NEG	POS
12	POS	NEG	POS	NEG	POS
13	NEG	NEG	NEG	NEG	NEG
14	POS	NEG	NEG	POS	POS
15	POS	NEG	POS	NEG	POS
16	NEG	NEG	NEG	POS	POS
17	NEG	NEG	NEG	POS	POS
18	POS	NEG	POS	NEG	POS
19	POS	NEG	POS	NEG	POS
20	NEG	NEG	NEG	NEG	POS
Test Date	15 Dec 08, 12 Jun 09		15 Dec 08, 10 Feb 09	13 Aug 09, 24 Aug 09	
Test Site	SC		SC	SC	
Kit Part Code	25228		39-6341	2647	
Kit Lot No.	956768, 960383		6611458, 6612575	090216	
Kit Exp. Date	21 Jul 09, 20 Feb 10		31 Mar 09, 30 Jun 09	22 Nov 09	
Kit Regulatory Status	IVD		IVD	CE	

¹Positive/reactive results are noted in bold red.

Panel members 5, 6, 9, 18 and 20 are no longer available. Data are provided for informational purposes only.

POS = Positive; NEG = Negative

IVD = In Vitro Diagnostic; CE = Conformité Européenne or CE Marking

AccuSet™ Anti-HIV-1 Low Titer Performance Panel Modified PRB109(M) (0800-0301) / Batch #10311491

HIV-1 Western Blot

Panel Member	Band Pattern	Calypse HIV-1 Western Blot	
		Result ^{1,2}	
01	24	IND	
02	No Bands	NEG	
03	f24	IND	
04	f24	IND	
05	f24, 160	IND	
06	f160	IND	
07	No Bands	NEG	
08	No Bands	NEG	
09	f24, f160	IND	
10	160	IND	
11	f24, f120, 160	IND	
12	24, f120, 160	POS	
13	18	IND	
14	18, f24, 160	IND	
15	18, 24, 160	POS	
16	f160	IND	
17	f24, 160	IND	
18	24, f41, f51, f120, 160	POS	
19	24, f41, f51, 55, 66, f120, 160	POS	
20	24, f55, 160	POS	
Test Date		29-May-09	
Test Site		SC	
Kit Part Code		98002	
Kit Lot No.		B9295-053	
Kit Exp. Date		22-Oct-09	
Kit Regulatory Status		IVD	

¹Western Blots were interpreted using ASTPHLD/CDC criteria (MMWR, vol. 38, S-7, 1989).

²Positive/reactive results are noted in bold red.

Panel members 5, 6, 9, 18 and 20 are no longer available. Data are provided for informational purposes only.

f = Faint; POS = Positive; IND= Indeterminate; NEG = Negative

IVD = In Vitro Diagnostic

The package insert for this panel can be found at www.seracare.com.

A printed copy of the package insert or data sheet may be requested by email at info@seracare.com, or by phone at 508.244.6400.