

AccuSet™ Toxoplasmosis Performance Panel

0820-0321 / Batch #10344475

OVERVIEW


AccuSet™ Toxoplasmosis Performance Panel (0820-0321 / Batch #10344475) is a 24-member validation panel of undiluted, naturally occurring plasma samples (1 vial per member, 1.0 mL per vial). Panel members represent bleeds from multiple individuals positive for antibodies to toxoplasmosis. Each sample represents a single collection event. No preservatives were added.

Test results from commercially-available toxoplasmosis assays are included for characterization of the panel members. This panel of human plasma samples demonstrates a range of antibody reactivity which approach the sensitivity limits for several toxoplasmosis test methods. One sample is included as a non-reactive sample and is negative for all toxoplasmosis test methods performed.

For Research Use Only. Not for use in diagnostic procedures. Data are offered for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.

CAUTION: The units that make up this panel were tested and found negative for anti-HIV-1/2, HBsAg and anti-HCV. This does not ensure the absence of these or other human pathogens.

AccuSet™ Toxoplasmosis Performance Panel



This graph demonstrates reactivity among panel members from the Trinity Biotech Captia™ Toxo IgM and Toxo IgG assays.

AccuSet™ Toxoplasmosis Performance Panel

0820-0321 / Batch #10344475

Panel Member Information

Panel Member	SeraCare Batch #	SeraCare Donor ID #	Bleed Date
01	10311431	BD360819	28-Sep-2015
02	10311432	BD360820	23-Oct-2015
03	10321177	BD357329	06-May-2016
04	10321178	BD357329	08-May-2016
05	10321179	BD357330	17-Jan-2018
06	10321180	BD357330	19-Jan-2018
07	10007625	BD200224	30-May-2013
08	9259176	BD110978	29-Apr-2012
09	10127362	BD250523	11-Jun-2015
10	10127410	BD250568	12-Jun-2015
11	9259496	BD106585	21-May-2012
12	10127394	BD250554	15-Jun-2015
13	10127482	BD250640	12-Jun-2015
14	10127353	BD250514	11-Jun-2015
15	10181204	BD281977	06-Jan-2016
16	10127380	BD250541	12-Jun-2015
17	10181175	BD281948	12-Feb-2016
18	10127387	BD250547	15-Jun-2015
19	10127483	BD250641	12-Jun-2015
20	10181186	BD281959	05-Feb-2016
21	10127501	BD250659	15-Jun-2015
22	10181161	BD281935	13-Feb-2016
23	10181176	BD282088	12-Feb-2016
24	10181216	BD281989	19-Feb-2016

AccuSet™ Toxoplasmosis Performance Panel

0820-0321 / Batch #10344475

Toxoplasmosis IgM

Panel Member	Trinity Biotech Captia™ Toxoplasma gondii IgM (s/co) ^{1,5}	Abbott ARCHITECT Toxo IgM (INDEX) ^{2,5}	bioMerieux VIDAS® Toxo IgM (TV) ^{3,5}	DiaSorin LIAISON® Toxo IgM (AU/mL) ^{4,5}	Bio-Rad Platelia™ Toxo IgM (s/co) ^{1,5}
01	1.7	2.6	1.9	12.9	2.8
02	1.2	1.9	2.0	16.1	3.1
03	1.8	1.4	2.8	75.7	3.3
04	1.6	1.2	2.8	68.1	3.2
05	2.6	4.0	2.7	71.1	3.3
06	3.4	2.8	2.3	45.7	2.9
07	3.4	1.0	1.0	11.5	1.7
08	0.7	0.8	0.3	<3.0	0.8
09	0.1	0.1	0.1	<3.0	0.0
10	0.6	0.1	0.2	6.7	0.2
11	0.1	0.1	0.1	<3.0	0.0
12	0.1	0.1	0.0	<3.0	0.0
13	0.2	0.3	0.5	<3.0	1.1
14	0.1	0.1	0.1	<3.0	0.0
15	0.4	0.4	0.5	3.1	0.8
16	0.1	0.2	0.3	<3.0	0.6
17	0.5	0.1	0.2	<3.0	0.2
18	0.4	0.1	0.1	<3.0	0.0
19	0.1	0.1	0.1	<3.0	0.0
20	0.3	0.0	0.0	<3.0	0.0
21	0.1	0.0	0.1	<3.0	0.0
22	0.2	0.1	0.1	3.3	0.0
23	0.7	0.1	0.2	<3.0, 3.5 ⁶	0.1
24	0.3	0.1	0.1	<3.0	0.0
Test Date	22-May-2018	28-May-2018	23-May-2018	23-May-2018	24-May-2018
Test Site	SC	RL	RL	RL	SC
Kit Part Code	2325160	6C20-25	NA	NA	26211
Kit Lot No.	2325160-085	84166L100	1006192440	NA	7F0035
Kit Exp. Date	31-Mar-2019	27-Nov-2018	09-Oct-2018	NA	15-Oct-2018
Kit Regulatory Status	IVD/CE	IVD/CE	IVD	IVD/CE	IVD/CE

¹Results are reported as a signal to cutoff ratio (s/co); positive/reactive results are noted in bold red.

²Results are reported as an index value; positive/reactive results are noted in bold red.

³Results are reported as a test value (TV); positive/reactive results are noted in bold red.

⁴Results are reported as an arbitrary unit (AU/mL); positive/reactive results are noted in bold red.

⁵Results are reported as the mean result of duplicate testing.

⁶Results are reported as individual replicates.

SC = SeraCare; RL = Reference Lab; NA = Not Available

IVD = In Vitro Diagnostic; CE = Conformité Européenne or CE Marking

AccuSet™ Toxoplasmosis Performance Panel

0820-0321 / Batch #10344475

Toxoplasmosis IgG

Panel Member	Abbott ARCHITECT Toxo IgG (IU/mL) ^{1,3}	bioMerieux VIDAS® Toxo IgG II (IU/mL) ^{1,3}	DiaSorin LIAISON® Toxo IgG (IU/mL) ^{1,3}	Trinity Biotech Captia™ Toxoplasma gondii IgG (s/co) ^{2,3}
01	19.3	135.5	70.2	5.1
02	65.2	281.0	125.0	5.1
03	122.0	263.5	207.5	4.9
04	104.4	269.5	190.0	4.7
05	24.3	43.0	32.8	2.9
06	24.9	45.5	34.1	3.2
07	171.7	>300	282.5	5.9
08	24.0	141.0	132.5	5.9
09	4.8	29.0	19.3	3.1
10	19.3	192.0	73.4	5.1
11	0.0	0.0	0.0	0.1
12	7.9	36.0	19.4	3.3
13	1.8	18.5	19.0	2.7
14	5.9	41.5	34.7	3.5
15	14.3	76.5	40.4	4.2
16	38.1	126.5	116.5	5.5
17	66.7	>300	200.0	5.7
18	20.3	176.5	98.8	5.1
19	7.6	40.0	30.0	3.3
20	1.2	16.5	10.6	2.7
21	6.3	39.0	35.6	3.2
22	8.9	66.5	36.2	4.0
23	21.1	88.0	46.6	4.2
24	42.8	>300	148.0	5.6
Test Date	28-May-2018	22-May-2018	23-May-2018	22-May-2018
Test Site	RL	RL	RL	SC
Kit Part Code	6C19-25	NA	NA	2325100
Kit Lot No.	84096L100	1006162530	NA	2325100-580
Kit Exp. Date	19-Dec-2018	28-Oct-2018	NA	31-May-2019
Kit Regulatory Status	IVD/CE	IVD	IVD	IVD/CE

¹Results are reported as international units (IU/mL); positive/reactive results are noted in bold red.

²Results are reported as a signal to cutoff ratio (s/co); positive/reactive results are noted in bold red.

³Results are reported as the mean result of duplicate testing.

SC = SeraCare; RL = Reference Lab; NA = Not Available
IVD = In Vitro Diagnostic; CE = Conformité Européenne or CE
Marking

The package insert for this panel can be found at
www.seracare.com.

A printed copy of the package insert or data sheet
may be requested by email at info@seracare.com or
by phone at 508.244.6400.