

AccuSpan™ HCV RNA Linearity Panel

PHW805 (2410-0166) / Batch #10408955

OVERVIEW


The AccuSpan™ HCV RNA Linearity Panel PHW805 (2410-0166) / Batch #10408955 is an eight-member panel made from serial dilutions of high titer HCV positive plasma with established reactivity for HCV (Hepatitis C) RNA. This panel consists of seven members representing serial log dilutions of HCV positive plasma in HCV RNA negative diluent and one negative member prepared from the diluent. The diluent was prepared from normal human plasma that was filtered through a 0.2 micron filter. Sodium azide (0.09%) was added as a preservative.

Results are reported for each panel member on specific test methods. Linearity is shown graphically by plotting observed results against expected results. A series of dilutions from the WHO International Standard was tested in the same run as the AccuSpan HCV RNA Linearity Panel members. Both expected and observed results from the dilutions series are reported; the expected values from the WHO standards are the dilution estimated values.

For Research Use Only. Not for use in diagnostic procedures. Data are provided for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.

CAUTION: These materials have not been treated and should be considered biohazardous. Follow Universal Precautions.

AccuSpan™ HCV RNA Linearity Panel



HCV RNA results were obtained using the Roche COBAS® AmpliPrep/COBAS® TaqMan® HCV Version 2.0 test method. Member #1 is the mean result of triplicate testing from a dilution of 1:10 and 1:100. Member #2 - #6 is the mean of three replicates. Member #7 is the average of 2 results as one replicate tested below the assay limit of quantitation. A line of best fit is shown.

AccuSpan™ HCV RNA
Linearity Panel
PHW805 (2410-0166) / Batch #10408955

HCV RNA

Panel Member	Roche COBAS® AmpliPrep/COBAS® TaqMan® HCV Version 2.0 ^{1,2} (log IU/mL)	Abbott m2000 RealTime HCV ^{1,2} (log IU/mL)
01	7.25 ³	7.20
02	6.17	6.25
03	5.28	5.30
04	4.37	4.28
05	3.46	3.26
06	2.50	2.35
07	1.42 ⁴ , <15	1.42
08	TND	TND
Test Date	02-May-2019	02-May-2019
Test Site	RL	RL
Test Kit Range	15 to 100,000,000 IU/mL 1.18 to 8.00 log IU/mL	12 to 100,000,000 IU/mL 1.08 to 8.00 log IU/mL
Test Kit Part Code	NA	NA
Test Kit Lot No.	E2156800000	NA

¹Results are reported as log International Units per mL (log IU/mL); positive/reactive results are noted in bold red.

²Results are reported as the mean result of triplicate testing.

³Results are reported as the mean result of triplicate testing from a dilution of 1:10 and 1:100.

⁴Results are reported as the mean result of two replicates, one replicate tested below the assays limit of quantitation and the result is displayed independently.

RL = Reference Lab; NA = Not available

AccuSpan™ HCV RNA
Linearity Panel
PHW805 (2410-0166) / Batch #10408955

5th WHO International HCV RNA Standard (14/150)

Sample ID	Expected Values (log IU/mL)	Observed Values on Roche COBAS® AmpliPrep/COBAS® TaqMan® HCV Version 2.0	
		(log IU/mL) ¹	% Difference ²
Sample 1	4.60	4.44	-3.58
Sample 2	4.00	3.92	-1.97
Sample 3	3.70	3.96	6.94
Sample 4	3.00	3.22	7.43
Test Date		02-May-2019	
Test Site		RL	
Test Kit Range		15 to 100,000,000 IU/mL 1.18 to 8.00 log IU/mL	
Kit Part Code		NA	
Kit Lot No.		E2156800000	

¹WHO panel was tested in the same test run as the AccuSpan™ HCV RNA Linearity Panel members. Samples were run in singlet. Positive/reactive results are noted in bold red.

²Percentage difference is how much the observed concentration differs from the expected concentration. Values calculated for reference only. Laboratories may use the data to apply a correction factor to the test results.

RL = Reference Lab; NA = Not available

The package insert for this panel can be found at www.seracare.com.

A printed copy of the package insert or data sheet may be requested by email at info@seracare.com or by phone at 508.244.6400.