

KPL Propidium Iodide

Item No.

5930-0008 (71-04-01)

Size

1.0 mg

DESCRIPTION

KPL Propidium Iodide (PI) is a fluorescent nucleic acid dye which binds only to double-stranded nucleic acids. As a nuclear counterstain, propidium iodide stains the nucleus light to dark red in color. The fluorescent red color provides excellent contrast to green or blue fluorochromes such as FITC, Cy2 or AMCA.

FORM/STORAGE

Propidium Iodide is supplied as a dark red powder. Store powder at 2-8°C. Stable for a minimum of 1 year from date of receipt when stored at 2-8°C.

CONTENT

1 x 1.0 mg of KPL Propidium Iodide 5930-0007 (71-04-00).

RECONSTITUTION

Reconstitute propidium iodide powder with 3.0 mL of PBS to a final concentration of 0.33 mg/mL. The reconstituted stock solution should be stored at -20°C and is stable for approximately 6 months.

APPLICATIONS

KPL Propidium Iodide can be used for fluorescent *in situ* hybridization, immunohistochemistry and flow cytometry.

SUGGESTED PROTOCOL

The dilutions and staining times stated here are offered as suggestions. Due to individual preference for stain intensity, these parameters should be optimized by the user.

Fluorescent *in Situ*

Hybridization/Immunohistochemistry:

1. Detect sample using an appropriate fluorescent conjugate.
2. Dilute the KPL Propidium Iodide stock solution 1:1000-1:30,000 in 1X PBS.
3. Apply diluted KPL Propidium Iodide generously to the slide and incubate at room temperature for 3 - 5 minutes in the dark (cover slide with aluminum foil).
4. Rinse the slide for 5 - 10 seconds in 1X PBS.
5. Mount with KPL Fluorescent Mounting Media or other appropriate mounting media.

NOTE: Intense KPL Propidium Iodide staining may mask other florescent signals. It is recommended that the concentration of KPL PI be titrated to determine the optimal dilution.

SOLUTION PREPARATION

PBS: Add 8.0 g NaCl, 0.2 g KCl, 1.44 g Na₂HPO₄ and 0.24 g KH₂PO₄ to 500 mL of reagent quality water. Adjust the pH to 7.4 with HCl. QS to 1 Liter with reagent quality water.

PRODUCT SAFETY AND HANDLING

See SDS (Safety Data Sheet) for this product.

RELATED PRODUCTS

	CAT. NO
KPL Fluorescent Mounting Media	5570-0005 (71-00-16)
KPL DAPI	5930-0006 (71-03-01)

Go to www.seracare.com for a complete list of fluorochrome labeled antibodies and streptavidin.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.