

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Klinefelter Syndrome (XXY) Reference Material

Material Number: 0720-0953 Batch Number: 10620439

Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).

Concentration
(Qubit dsDNA BR Assay): 28.6 ng/mL

Fill Volume: 1.0 mL

Date of Manufacture: 06-MAY-2022 Expiration Date: 04-MAR-2026


Method used for NIPT: Sage NIPT Test
LDT version of Veriseq solution v2 (illumina Inc)

Fetal Fraction: 8.7 %

Call for XXY: Pass

Storage: 2-8 °C

Approval:



12 MAY 2022

Prepared By

Date