

Cell-Ref™ FFPE Human Cancer Cell Slide - CCRF-CEM

Catalog Number:	ASO-1056
Size:	1 pack (5 slides)
Format:	FFPE
Volume/Dimension:	5 µm thick section with 7.0mm diameter cell core
Quantity:	200,000 cells
Concentration:	200,000 cells/slide
Cell(s):	CCRF-CEM
Tissue Type:	Leukemia
Shelf Life:	24 months
Shipping Temperature:	Ambient
Fixation Process:	10% buffered formalin solution (4% formaldehyde); 24 hrs; room temperature
Storage Temperature:	-20°C
Quality Control:	Quality: > 60% Cell Density (H&E Staining of 5µm cell slide)
Regulatory:	Research Use Only. Not Intended for human or animal diagnostic or therapeutic use.
Restrictions:	This is a limited use product. This product, any material that contains this product in whole or in part, any progeny, modification or derivative of this product, any cell or animal made or modified by using this product, and any progeny, modification or derivative of such cell or animal may not be transferred by the purchaser to any other person, entity or any of the above to perform services for the benefit of any other person or entity or for commercial purpose of any kind. This product may only be used by the purchaser for its internal research for use as a research tool for research purposes.

DS-ASO-1056-V2