

## Cell-Ref™ FFPE Human Cancer Cell Slide - SUDHL6

|                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number:</b>       | <b>ASO-1060</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Size:</b>                 | 1 pack (5 slides)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Format:</b>               | FFPE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Volume/Dimension:</b>     | 5 µm thick section with 7.0mm diameter cell core                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Quantity:</b>             | 200,000 cells                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Concentration:</b>        | 200,000 cells/slide                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Cell(s):</b>              | SUDHL6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Tissue Type:</b>          | Lymphoma                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Shelf Life:</b>           | 24 months                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Shipping Temperature:</b> | Ambient                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Fixation Process:</b>     | 10% buffered formalin solution (4% formaldehyde); 24 hrs; room temperature                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage Temperature:</b>  | -20°C                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Quality Control:</b>      | Quality: > 60% Cell Density (H&E Staining of 5µm cell slide)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Regulatory:</b>           | Research Use Only. Not Intended for human or animal diagnostic or therapeutic use.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Restrictions:</b>         | This is a limited use product. This product, any material that contains this product in whole or in part, any progeny, modification or derivative of this product, any cell or animal made or modified by using this product, and any progeny, modification or derivative of such cell or animal may not be transferred by the purchaser to any other person, entity or any of the above to perform services for the benefit of any other person or entity or for commercial purpose of any kind. This product may only be used by the purchaser for its internal research for use as a research tool for research purposes. |

DS-ASO-1060-V2