

## Cell-Ref™ FFPE Human Cancer Cell Slide - CHO-K1

<b>Catalog Number:</b>	<b>ASO-1072</b>
<b>Size:</b>	1 pack (5 slides)
<b>Format:</b>	FFPE
<b>Volume/Dimension:</b>	5 µm thick section with 7.0mm diameter cell core
<b>Quantity:</b>	200,000 cells
<b>Concentration:</b>	200,000 cells/slide
<b>Cell(s):</b>	CHO-K1 (hamster)
<b>Tissue Type:</b>	Hamster
<b>Shelf Life:</b>	24 months
<b>Shipping Temperature:</b>	Ambient
<b>Fixation Process:</b>	10% buffered formalin solution (4% formaldehyde); 24 hrs; room temperature
<b>Storage Temperature:</b>	-20°C
<b>Quality Control:</b>	Quality: > 60% Cell Density (H&E Staining of 5µm cell slide)
<b>Regulatory:</b>	For Research Use Only. Not intended for human or animal diagnostics or therapeutics.
<b>Restrictions:</b>	This is a limited use product. This product, any material that contains this product in whole or in part, any progeny, modification or derivative of this product, any cell or animal made or modified by using this product, and any progeny, modification or derivative of such cell or animal may not be transferred by the purchaser to any other person, entity or any of the above to perform services for the benefit of any other person or entity or for commercial purpose of any kind. This product may only be used by the purchaser for its internal research for use as a research tool for research purposes.

DS-ASO-1072-V2