

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq Breast CNV Mix, +3 copies

Material No: 0710-0411 Batch Number: 10311352

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs


Date of Manufacture: 30Jan2018 Expiration Date: 30Jan2020

Test Method: Qubit dsDNA BR Assay, BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System.

Concentration: Nominal value: 10 ng/μL
 Average measured value: 12.30 ng/μL
 Fill volume: 20 μL

| Gene | Number of Additional Copies by Digital PCR |
|-------|--|
| ERBB2 | 3.68 |
| FGFR3 | 3.65 |
| MYC | 3.77 |

Approval:



Prepared By

02/05/2018
 Date