

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 21 Male - Matched Reference Material

Material No: 0720-0167

Batch Number: 10417417

Material Description: Mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).

Date of Manufacture: 10-May-2019

Expiration Date: 17-Mar-2023

Concentration (Qubit dsDNA
BR Assay): 30.9 ng/mL

Volume: 1.0 mL

Fetal Fraction by NGS
(Illumina VeriSeq): 11%

Trisomy Call for
Chromosome 21 by NGS
(Illumina VeriSeq): Pass

Storage: 2-8 °C

Approval:

A handwritten signature in black ink, appearing to read "J. [unclear]".

Prepared By

21-May-2019

Date