


# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® Euploid Female - Matched Reference Material  
Material Number: 0720-0170 Batch Number: 10494998  
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).  
Concentration (Qubit dsDNA BR Assay): 34.0 ng/mL  
Volume: 1.0 mL  
Date of Manufacture: 11 MAY 2020 Expiration Date: 11 MAY 2024  
Fetal Fraction by NGS (Innatal® Prenatal Screen): 11.8%  
  
Euploid Call by NGS (Innatal Prenatal Screen): Pass  
Call for fetal gender by NGS (Innatal Prenatal Screen): XX  
Storage: 2-8 °C

**Approval:**

A handwritten signature in blue ink, appearing to be "R. A. ...".

Prepared By

26 May 2020

Date