

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® Euploid Female - Matched Reference Material  
Material Number: 0720-0170 Batch Number: 10643581  
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).  
Concentration (Qubit dsDNA BR Assay): 31.5 ng/mL  
Volume: 1.0 mL  
Date of Manufacture: 21 NOV 2022 Expiration Date: 21 NOV 2026  
Fetal Fraction by NGS (Illumina VeriSeq): 18.0%  
  
Euploid Call by NGS (Illumina VeriSeq): Pass  
Call for fetal gender by NGS (Illumina VeriSeq): XX  
Storage: 2-8 °C

**Approval:**

A handwritten signature in blue ink, appearing to be "MA".

05-DEC-2022

Prepared By

Date