

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Mutation Mix AF0.5%

Material Number: 0710-0531 Batch Number: 10646979

Material Description: A ctDNA-like mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Concentration
(Qubit dsDNA BR Assay): Nominal value: 10 ng/μL; Average measured value: 15.6 ng/μL

Fill Volume: 25 μL

Date of Manufacture: 13 DEC 2022 Expiration Date: 13 DEC 2024

Storage: -20°C

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Mutation Mix AF0.5%

Digital PCR testing
using BioRad QX200™
Droplet Digital™ PCR
System:

Gene ID	COSMIC Identifier	Amino Acid Change	Average AF%
AKT1	COSM33765	p.E17K	0.51
BRAF	COSM476	p.V600E	0.48
EGFR	COSM6224	p.L858R	0.54
EGFR	COSM6240	p.T790M	0.46
ERBB2	COSM20959	p.A775_G776insYVMA	0.42
KIT	COSM1314	p.D816V	0.53
KRAS	COSM521	p.G12D	0.56
NCOA4/RET	NA	Translocation	0.56
NRAS	COSM584	p.Q61R	0.54
PIK3CA	COSM775	p.H1047R	0.55
PIK3CA	COSM12464 ¹	p.N1068fs*4	0.55
EML4-ALK	NA	Translocation	0.49
ALK	COSM144250	p.G1202R	0.56
ALK	COSM28055	p.F1174L	0.56
BRCA1	COSM1383519	p.K654fs*47	0.46
BRCA2	COSM1738242	p.R2645fs*3	0.50
EGFR	COSM12370	p.L747_P753>S	0.59
EGFR	COSM6256	p.S752_I759delSPKANKEI	0.50
EGFR	COSM6223	p.E746_A750delELREA	0.64
KRAS	COSM516	p.G12C	0.52
CD74/ROS1	NA	Translocation	0.53
KRAS	COSM554	p.Q61H	0.56

Gene ID	Average CNV in ctDNA ²	Average Additional Copies (per cell) in ctDNA
ERBB2	2.56	0.56
MET	2.41	0.41
MYC	2.37	0.37

NA = not applicable

¹As of June 2019, this mutation is no longer listed in the COSMIC database.

²Compare to a normal CN of 2.00.

Next Generation Sequencing testing using Archer® Reveal ctDNA™ 28 Kit run on an Illumina® MiSeq™ using v2 (2x150 bp) PE chemistry reagents¹:

Gene ID	COSMIC Identifier	Amino Acid Change	AF%
AKT1	COSM33765	p.E17K	0.54
BRAF	COSM476	p.V600E	0.54
EGFR	COSM6224	p.L858R	0.50
EGFR	COSM6240	p.T790M	0.46
ERBB2	COSM20959	p.A775_G776insYVMA	0.31
KIT	COSM1314	p.D816V	0.67
KRAS	COSM521	p.G12D	0.70
NCOA4/RET	NA	Translocation	NA
NRAS	COSM584	p.Q61R	0.72
PIK3CA	COSM775	p.H1047R	0.59
PIK3CA	COSM12464 ³	p.N1068fs*4	0.74
EML4-ALK	NA	Translocation	NA
ALK	COSM144250	p.G1202R	0.40
ALK	COSM28055	p.F1174L	0.58
BRCA1	COSM1383519	p.K654fs*47	NA
BRCA2	COSM1738242	p.R2645fs*3	NA
EGFR	COSM12370	p.L747_P753>S	0.51
EGFR	COSM6256	p.S752_I759delSPKANKEI	0.47
EGFR	COSM6223	p.E746_A750delELREA	0.65
KRAS	COSM516	p.G12C	0.70
CD74/ROS1	NA	Translocation	NA
KRAS	COSM554	p.Q61H	0.42

Gene ID	CNV in ctDNA ⁴	Additional Copies (per cell) in ctDNA
ERBB2	1.29	0.58
MET	1.27	0.54
MYC	NA	NA

NA = not applicable; AF% and CNV marked NA were not targeted by the panel.

¹NGS was performed as an orthogonal verification step. Parameters used:

DNA input = 50 ng

of samples / flow cell = 2

of total reads / sample = 4.7M

Average read depth = 8548

On-target reads = 96.4%

Q30 score = 89.3%

Analysis = Archer Analysis Suite v6.2.7 (default settings)

²Please see the poster from NIST for more information about assay sensitivity:

<https://digital.seracare.com/multilab-assessment-reference-materials-ctdna-poster2018>

³As of June 2019, this mutation is no longer listed in the COSMIC database.

⁴Compare to a normal CN of 2.00.

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Mutation Mix AF0.5%

Note: The MET gene is amplified using two synthetic constructs with a small region of overlap between the constructs (see package insert for genomic coordinates). Assays which target this region of overlap may report higher amplification levels.

Approval:

A handwritten signature in black ink, appearing to be "J. M. ...".

15 DEC 2022

Prepared By

Date