

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Reference Material AF 5%

Material Number: 0710-0669 Batch Number: 10659790

Material Description: A ctDNA-like mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Concentration (Qubit dsDNA BR Assay): Nominal value: 25 ng/mL; Average measured value after extraction using Qiagen QIAamp Circulating Nucleic Acid Kit: 32.2 ng/mL

Fill Volume: 5.0 mL

Date of Manufacture: 05 APR 2023 Expiration Date: 05 APR 2027

Storage: 2-8°C

Digital PCR testing using
BioRad QX200™ Droplet
Digital™ PCR System:

Gene ID	COSMIC Identifier	Amino Acid Change	Average AF%
AKT1	COSM33765	p.E17K	5.22
BRAF	COSM476	p.V600E	5.26
EGFR	COSM6224	p.L858R	5.30
EGFR	COSM6240	p.T790M	5.03
ERBB2	COSM20959	p.A775_G776insYVMA	5.12
KIT	COSM1314	p.D816V	5.18
KRAS	COSM521	p.G12D	5.35
NCOA4/RET	NA	Translocation	5.07
NRAS	COSM584	p.Q61R	5.47
PIK3CA	COSM775	p.H1047R	4.98
PIK3CA	COSM12464 ¹	p.N1068fs*4	4.98
EML4-ALK	NA	Translocation	4.58
ALK	COSM144250	p.G1202R	4.88
ALK	COSM28055	p.F1174L	4.88
BRCA1	COSM1383519	p.K654fs*47	4.79
BRCA2	COSM1738242	p.R2645fs*3	4.61
EGFR	COSM12370	p.L747_P753>S	6.02
EGFR	COSM6256	p.S752_I759delSPKANKEI	5.14
EGFR	COSM6223	p.E746_A750delELREA	5.72
KRAS	COSM516	p.G12C	5.28
CD74/ROS1	NA	Translocation	4.88
KRAS	COSM554	p.Q61H	4.89

Gene ID	Average CNV in ctDNA ²	Average Additional Copies (per cell) in ctDNA
ERBB2	7.67	5.67
MET	6.11	4.11
MYC	6.27	4.27

NA = not applicable

¹As of June 2019, this mutation is no longer listed in the COSMIC database.

²Compare to a normal CNV of 2.00.

Next Generation Sequencing testing using Archer® Reveal ctDNA™ 28 Kit run on an Illumina® MiSeq™ using v2 (2x150 bp) PE chemistry reagents^{1,2}:

Gene ID	COSMIC Identifier	Amino Acid Change	AF%
AKT1	COSM33765	p.E17K	5.00
BRAF	COSM476	p.V600E	5.00
EGFR	COSM6224	p.L858R	4.00
EGFR	COSM6240	p.T790M	4.80
ERBB2	COSM20959	p.A775_G776insYVMA	4.10
KIT	COSM1314	p.D816V	5.90
KRAS	COSM521	p.G12D	5.90
NCOA4/RET	NA	Translocation	NA
NRAS	COSM584	p.Q61R	5.70
PIK3CA	COSM775	p.H1047R	5.30
PIK3CA	COSM12464 ³	p.N1068fs*4	5.20
EML4-ALK	NA	Translocation	NA
ALK	COSM144250	p.G1202R	4.70
ALK	COSM28055	p.F1174L	4.70
BRCA1	COSM1383519	p.K654fs*47	NA
BRCA2	COSM1738242	p.R2645fs*3	NA
EGFR	COSM12370	p.L747_P753>S	5.20
EGFR	COSM6256	p.S752_I759delSPKANKEI	4.50
EGFR	COSM6223	p.E746_A750delELREA	4.50
KRAS	COSM516	p.G12C	5.00
CD74/ROS1	NA	Translocation	NA
KRAS	COSM554	p.Q61H	4.30

Gene ID	CNV in ctDNA ⁴	Additional Copies (per cell) in ctDNA
ERBB2	3.11	4.22
MET	3.56	5.12
MYC	NA	NA

NA = not applicable; AF% and CNV marked NA were not targeted by the panel.

¹NGS was performed as an orthogonal verification step. Parameters used:

DNA input = 50 ng

of samples / flow cell = 7

of total reads / sample = 5.5 M

Average read depth = 10368X

On-target reads = 96.7%

Q30 score = 91.49%

Analysis = Archer Analysis Suite v6.2.7 (default settings except for: Error correction)

²Please see the poster from NIST for more information about assay sensitivity:

<https://digital.seracare.com/multilab-assessment-reference-materials-ctdna-poster2018>

³As of June 2019, this mutation is no longer listed in the COSMIC database.

⁴Compare to a normal CNV of 2.00.

Note: The MET gene is amplified using two synthetic constructs with a small region of overlap between the constructs (see package insert for genomic coordinates). Assays which target this region of overlap may report higher amplification levels.

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Reference Material AF 5%

Approval:

A handwritten signature in blue ink, appearing to be "M A".

26 APR 2023

Prepared By

Date