

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® TNA WT Mix

Material Number: 0710-1580 Batch Number: 10500560

Material Description: Seraseq® TNA WT Mix contains DNA (0710-1582) and RNA (0710-1584) from the reference human cell line GM24385. The cell line is WT in that it is not derived from cancerous tissue, but mutations from the human reference sequence (GrCh38) are still present.

Fill Volume: 20 µL

Date of Manufacture: 18 JUN 2020 Expiration Date: 16 JUN 2022


Test Method for Concentration: Qubit dsDNA BR Assay and Qubit RNA BR Assay

Concentration of DNA Vial: 28.7 ng/µL Concentration of RNA Vial: 31.1 ng/µL

Test Method for Variant Detection: Illumina TruSight Oncology 500 Assay

Variant Detection Results: No pathogenic variants identified. VCF files available upon request. TMB Score of 0 reported.

Approval:



Prepared By

19 Jun 2020
Date