

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® Trisomy 21 Female - Matched Reference Material

Material Number: 0720-0168      Batch Number: 10456868

Material Description: A ccfdNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).

Concentration  
(Qubit dsDNA BR Assay): 31.2 ng/mL

Volume: 1.0 mL

Date of Manufacture: 11 Dec 2019      Expiration Date: 19 Nov 2023

Fetal Fraction by NGS  
(Illumina VeriSeq): 13%

Trisomy Call for  
Chromosome 21 by NGS  
(Illumina VeriSeq): Pass

Call for fetal gender by  
NGS (Illumina VeriSeq): XX

Storage: 2-8°C

Approval:

A handwritten signature in black ink, appearing to be "J. Blum", written over a horizontal line.

Prepared By

16 DEC 2019

Date