

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® Trisomy 21 Female - Matched Reference Material  
 Material Number: 0720-0168 Batch Number: 10571861  
 Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).  
 Concentration (Qubit dsDNA BR Assay): 34 ng/mL  
 Volume: 1.0 mL  
 Date of Manufacture: 06 AUG 2021 Expiration Date: 06 JUL 2025  
 Fetal Fraction by NGS (Illumina VeriSeq):  
 Enter NGS data and type of NGS assay in gray. More lines may be added depending on how many different types of NGS methods were used. Some data, like fetal gender, may not be available from NGS data. If that is the case, delete those lines. Delete THIS  
 8.1%  
 Trisomy Call for Chromosome 21 by NGS (Illumina VeriSeq): Pass  
 Call for fetal gender by NGS (Illumina VeriSeq): XX  
 Storage: 2-8°C

Approval:



Prepared By

10 AUG 2021  
Date