

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq<sup>®</sup> PGT-A Trisomy 13 Reference Material

Material Number: 0720-0777                      Batch Number: 10683804

Material Description: Genomic DNA extracted from a male human cell line carrying three copies of chromosome 13.

Concentration (Qubit dsDNA HS Assay): 27.4 pg/μL

Fill Volume: 10 μL

Date of Manufacture: 23 OCT 2023                      Expiration Date: 25 AUG 2024

Trisomy Call using commercially available NGS-based PGT-A: Trisomy 13

Sex determined using commercially available NGS-based PGT-A: Male

Storage: -20 °C

**Approval:**

A handwritten signature in black ink, appearing to be "M. H. A.".

26 OCT 2023

Prepared By

Date